

To: «EDITOR2» «EDITOR3»

Publication: «PUBLICATION»



NEW PRODUCT RELEASE

NEW Dual Rate Syringe Pump

KD Scientific has introduced the **NEW Gemini 88 Syringe Pump** which is actually "Two Pumps in One".

The **Gemini 88** has two independent pumping channels linked through hardware and software.

The pump can infuse simultaneously at different rates, or infuse with one syringe and withdraw with the other. When combined with a valve box the **Model 88** can provide continuous delivery.

The **Gemini 88** is designed for continuous operation 24/7 with a wide range of syringes from 0.5 µl up to 140 ml.

The pump has high pressure capabilities and TTL and RS-232 interfaces for data acquisition and control. The communications ports enable daisy-chaining up to 100 pumps.

For complete specifications for the Gemini 88 contact KD Scientific or go to the KD Scientific website at www.kdscientific.com.

KD Scientific designs, manufactures and sells a range of quality fluidics equipment used by research laboratory markets worldwide.

KD Scientific syringe pumps are an economical solution to delivering precise and smooth flow in research, pilot plants and production applications. They are recognized worldwide for quality, accuracy and reliability. A broad line of syringe pumps are offered: from a simple one syringe infuse only, to a programmable multi-syringe infuse/withdrawal pump. **KD Scientific** operates as a separate business within Harvard Bioscience at their Holliston, Massachusetts facility.

For more information contact:

KD Scientific

Click on Image for 300 dpi jpg



84 October Hill Road
Holliston, Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: info@kdscientific.com

Web: www.kdscientific.com

Product and company names listed are trademarks or trade names of their respective companies.

NOTE: This Press Release has been sent electronically. If it has been sent to the wrong individual, or if you no longer wish to receive KD Scientific Press Releases, please notify us at: PressRelease@optonline.net