

**K**D Scientific offers a selection of Hamilton® microliter and specialty syringes. Please visit our website or contact our technical support department for the latest offerings.

# State-of-the-Art GASTIGHT® Syringes

## Premium Line of GASTIGHT® Syringes

Hamilton syringes are considered the industry standard for precision fluid delivery. KD Scientific offers a selection of these precision syringes which, when combined with our syringe pumps, offer fluid delivery with unparalleled accuracy and precision.

The syringes are designed with reinforced syringe plungers, syringe barrels, and are GASTIGHT® syringes with syringe volumes from 50 µl to 100 ml.

All Hamilton™ glass syringes are autoclavable when disassembled.



### 1000 Series GASTIGHT® Syringes

Item No.	Description
72-1831	1 ml PTFE LUER LOCK
72-1832	2.5 ml PTFE LUER LOCK
72-1833	5 ml PTFE LUER LOCK
72-1834	10 ml PTFE LUER LOCK
72-1835	25 ml PTFE LUER LOCK
72-1836	50 ml PTFE LUER LOCK
72-1837	100 ml PTFE LUER LOCK

### 1700 Series GASTIGHT® Syringes

Item No.	Description
72-1781	50 µl PTFE LUER LOCK
72-1782	100 µl PTFE LUER LOCK
72-1783	250 µl PTFE LUER LOCK
72-1784	500 µl PTFE LUER LOCK

GASTIGHT® is a registered trademark of the Hamilton Company.