

To: «EDITOR2» «EDITOR3»
Publication: «PUBLICATION»



NEW PRODUCT RELEASE

NEW Economical Syringe Pump with Smooth Flow!

KD Scientific has introduced the KDS 210 Syringe Pump which offers smooth flow with more advanced features than any other infusion/withdrawal pump in its price range.

There are five operating modes available consisting of Infusion, Withdrawal, Infusion then Withdrawal and Withdrawal then Infusion.

The **KDS 210** has independent rate and volume settings for infusion and withdrawal plus independent rate and volume settings for both infusion and withdrawal.

The **KDS 210** holds two syringes, 10 µl to 140 ml each. The units have built-in RS232C interface for computer "daisy chaining" up to 100 pumps.



KD Scientific designs, manufactures and sells a range of quality fluidics equipment used by research laboratory markets worldwide.

For more information contact:

KD Scientific

84 October Hill Road
Holliston , Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: info@kdsscientific.com

Web: www.kdsscientific.com

NOTE: This Press Release has been sent electronically. If it has been sent to the wrong individual, or if you no longer wish to receive KD Scientific Press Releases, please notify us at: PressRelease@optonline.net