

To: «EDITOR2» «EDITOR3»
Publication: «PUBLICATION»



NEW PRODUCT RELEASE

Ultra Smooth Flow Pump Perfect for MS Calibration

The KD Scientific Model KDS 100 is a basic, single syringe, infusion pump which combines precision with simplicity.

Any single syringe from 10 μ l to 60 ml is all that is needed for ESI sample introduction as well as calibration solutions.

The **KDS 100** is easy-to-use and has a simple menu-driven setup: select syringe size from displayed table, set dispense volume, and set dispense flow rate.

The **KDS 100** can do simple infusions or deliver a set volume. In the volume dispense mode, the **KDS 100** will deliver a specified volume and automatically shut off when that volume is delivered. When using this mode the volume delivered is continuously displayed on the LED until the delivery is complete. The advanced drive mechanism assures the smoothest most accurate flow for all delivery modes.

KD Scientific designs, manufactures and sells a wide range of quality fluidics equipment used by research laboratory markets worldwide.

For more information contact:

KD Scientific

84 October Hill Road
Holliston, Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: info@kdscientific.com

Web: www.kdscientific.com

Product and company names listed are trademarks or trade names of their respective companies.

NOTE: This Press Release has been sent electronically. If it has been sent to the wrong individual, or if you no longer wish to receive KD Scientific Press Releases, please notify us at: PressRelease@optonline.net



