



NEW PRODUCT RELEASE

New LabVIEW Driver for KD Scientific Products

KD Scientific products now can be powered by **National Instruments LabVIEW software**, a widely used graphical development environment for performing data acquisition analysis and presentation.

KD Scientific and National Instruments collaborated to develop the driver. This **NEW LabVIEW driver** enhances the automation and integration of the **KD 200 Series of Syringe Pumps** with other lab and process equipment in a variety of R&D and process applications.

The driver is a part of more than 4,900 new drivers from National Instruments that deliver connectivity for instruments serving many industries and all standard software and bus interfaces.

This **NEW Driver** demonstrates a commitment to connectivity and automation of lab equipment. Because engineers and scientists can view and modify the source code of these instrument drivers, they can easily integrate additional functionality and optimize code for their particular applications.

Engineers and scientists now can download **FREE** drivers for the **KD Scientific Series 200 Syringe Pumps** through the National Instrument's "**Instrument Driver Network**" at www.ni.com/idnet.

The Series 200 syringe pumps are widely used for delivering smooth consistent flow from 1 nl/hr to 145 ml/min. The units work with syringes from 10 ul to 140 ml. Typical accuracy is +/- 1%. There are 6 models in the Series 200 including the standard infusion only to a continuous pump.

Applications for the KDS 200 Series and these New Drivers are in pharmaceutical, chemical, petrochemical, biotechnology, semiconductor, plastics, industrial, government, and scientific research and development markets.

KD Scientific Inc. designs, manufactures and sells a wide range of quality fluidics equipment used by research laboratory markets worldwide.

For more information contact:

KD Scientific Inc.
84 October Hill Road



Holliston , Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: info@kdscientific.com

Web: www.kdscientific.com

Product and company names listed are trademarks or trade names of their respective companies.